



Since 2007

Master of Sterilization

HOSPICARE SALES & SERVICES

Deals in :- Hospital O.T. & CSSD Equipments

Fully Automatic Steam + ETO (COMBO) Sterilizer / Autoclave

Sailent Features

Steam Sterilizer

- Heating : Electrically Heated by Water Heater to Work 230 V 2000 Watt.
- Working Pressure : The Pressure Inside the Chamber is variable from 5 PSI to 20 PSI.
- Working Temp: 121°C Temperature is Fully controlled by PLC.
- Sterilization time 30-45 minutes.
- Fully Automatic Operated.

ETO Sterilizer

- Heating: Electrically heated by dry Heater to Work 230v & 400 Watts.
- Working Pressure: The Pressure inside the Chamber is variable from 5 PSI to 20 PSI
- Working Temp: 45°C to 50°C Constant.
- Suitable For 10-90, 20-80 EO-C02 Concentration.
- Perforated basket for ETO Materials.
- Sterilization time 4-5 Hours.
- Fully Automatic Operated.
- Vertical and Horizontal models available.



Common Sizes

16" x 24"

20" x 36"

20" x 48"

Manufacturing, Marketing & Services

HOSPICARE SALES & SERVICES

696, Gokul, Behind More Petrol Pump, Shivane, Pune - 411 023

Mobile No.:

9921101718 / 7020444869

Website:

www.hospicaresalesandservices.com

Email:

hospicarenitin01@gmail.com

sales@hospicaresalesandservices.com

service@hospicaresalesandservices.com

A Company That Cares For Innovation, Quality Products & Efficient After Sales Services