



Since 2007

Master of Sterilization

HOSPICARE SALES & SERVICES

Deals in :- Hospital O.T. & CSSD Equipments

Horizontal Steam Sterilizer / Autoclave

Sailent Features



- Chamber, Jacket, Door Locking Ring & Back Plate is made from S.S. 304
- Low water Level protection for heater,
- Two spring safety valves to prevent uncontrolled pressure.
- All measuring components.
- Vacuum system incorporated.
- Fulfill all requirement.
- Power consumption 440 V, 3-Phase.
- Best suited for rapid sterilization of instruments for Cardiac, Ortho, Neuro & Urology Of 20-50 Bedded Hospitals.
- Fully Automatic Operated, Semi Automatic Operated, & Manually Optional.
- Single Door & Double Door Options Available.

| | | | |
|------------------------|------------------|------------------|------------------|
| Common Sizes | 16" x 24" | 20" x 36" | 20" x 48" |
| Capacity in LTR | 75 LTR | 170 LTR | 280 LTR |

Manufacturing, Marketing & Services

HOSPICARE SALES & SERVICES

696, Gokul, Behind More Petrol Pump, Shivane, Pune - 411 023

Mobile No.:

9921101718 / 7020444869

Website:

www.hospicaresalesandservices.com

Email:

hospicarenitin01@gmail.com

sales@hospicaresalesandservices.com

service@hospicaresalesandservices.com

A Company That Cares For Innovation, Quality Products & Efficient After Sales Services