



Since 2007

Master of Sterilization

HOSPICARE SALES & SERVICES

Deals in :- Hospital O.T. & CSSD Equipments

Rectangular ETO Sterilizer / Autoclave



Sailent Features

- Chamber & outer cover made of Stainless Steel 304 Quality.
- Manual safety lock is provided on the door.
- Temperature indicator to indicate 0 Degree to 100 Degree.
- Pressure indicator to indicate plus and minus pressure.
- Incorporated vacuum pump to achieve vacuum in Chamber.
- Sterilization time 6-7 Hours.
- Fully Automatic Operated, Semi Automatic Operated, Manually Operated options available.

Common Sizes

12" x 12" x 24"

12" x 12" x 36"

12" x 12" x 48"

Manufacturing, Marketing & Services

HOSPICARE SALES & SERVICES

696, Gokul, Behind More Petrol Pump, Shivane, Pune - 411 023

Mobile No.:

9921101718 / 7020444869

Website:

www.hospicaresalesandservices.com

Email:

hospicarenitin01@gmail.com

sales@hospicaresalesandservices.com

service@hospicaresalesandservices.com

A Company That Cares For Innovation, Quality Products & Efficient After Sales Services