



Since 2007

Master of Sterilization

HOSPICARE SALES & SERVICES

Deals in :- Hospital O.T. & CSSD Equipments

Rectangular Steam Sterilizer / Autoclave



Sailent Features

- Best suited for parenteral solutions, injectables, hospital dressing And Linen.
- Operating Temperature 121°C or 134°C.
- All necessary safety precaution namely pressure locking Inbuilt.
- Incorporated Automatic Door Locking System.
- Loading carriage and trolley as standard accessories for easy transfer of load If Required.
- Choice of Single Door, Double Door, Flush Mounted, Panel Models.
- Steam Trap is incorporated in discharge line for fast and Perfect elimination of steam.
- Various material combination like Stainless Steel 304/316 is available.
- Fully Automatic Operated, Semi Automatic Operated, & Manually Operated Options Available.
- Best suited for more than 50 bedded hospitals.
- S.S.Flush Mounted Optional.

Standard Sizes	Capicity in LTR	Power Requirement	
600mm X 600mm X 1200mm Long	450 LTR	3 Phase	18 KW
600mm X 900mm X 1500mm Long	675 LTR	3 Phase	18 KW / 27 KW

Manufacturing, Marketing & Services

HOSPICARE SALES & SERVICES

696, Gokul, Behind More Petrol Pump, Shivane, Pune - 411 023

Mobile No.:

9921101718 / 7020444869

Website:

www.hospicaresalesandservices.com

Email:

hospicarenitin01@gmail.com

sales@hospicaresalesandservices.com

service@hospicaresalesandservices.com

A Company That Cares For Innovation, Quality Products & Efficient After Sales Services