



Since 2007

Master of Sterilization

# HOSPICARE SALES & SERVICES

Deals in :- Hospital O.T. & CSSD Equipments

## Vertical Steam Sterilizer / Autoclave



### Sailent Features

- Chamber & outer cover made of Stainless Steel 304 Quality.
- Lid is made of S.S 304 Quality & also silicon Gasket for its better life.
- 2 Nos. of Safety Valves.
- Compound Gauge is Incorporated to show pressure.
- Automatic Cut off Device is to protect the heater from running dry.
- Manual Operated & Automatic Operated Optional.

**Common Sizes**

**12" x 22"**

**16" x 24"**

Manufacturing, Marketing & Services By

## HOSPICARE SALES & SERVICES

696, Gokul, Behind More Petrol Pump, Shivane, Pune - 411 023

**Mobile No.:**

9921101718 / 7020444869

**Website:**

[www.hospicaresalesandservices.com](http://www.hospicaresalesandservices.com)

**Email:**

[hospicarenitin01@gmail.com](mailto:hospicarenitin01@gmail.com)

[sales@hospicaresalesandservices.com](mailto:sales@hospicaresalesandservices.com)

[service@hospicaresalesandservices.com](mailto:service@hospicaresalesandservices.com)

**A Company That Cares For Innovation, Quality Products & Efficient After Sales Services**